


[Kd Validated] Copb2 Rabbit Polymab® - P35606

Item Number: CM0004973



Introduction

Rabbit monoclonal antibody against human COPB2 (beta'-COP), validated for WB, IF/ICC, IHC-P, and ELISA. Detects human and mouse COPB2 (~102 kDa), a coatamer subunit essential for Golgi-to-ER retrograde transport.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] COPB2 Rabbit PolymAb®
Remarks/Alias	OPDD; MCPH19; beta'-COP
Species	Human
GeneID (Human)	9276
GeneID	9276
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 657-906 of human COPB2 (NP_004757.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P35606
Protein Weight	102kDa
Shipping	Ice bag