

Kd Validated Dnmt1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P26358

Item Number: CM0023649



Introduction

Rabbit polyclonal antibody to DNMT1 (DNA methyltransferase 1), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for DNA methylation and epigenetic inheritance studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] DNMT1 Rabbit pAb
Remarks/Alias	AIM; DNMT; MCMT; CXXC9; HSN1E; ADCADN; m.Hsal; [KD Validated] DNMT1
Species	Human
Gene ID (Human)	1786
Gene ID	1786
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-124 of DNMT1 (NP_001370.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26358
Protein Weight	183kDa
Shipping	Ice bag