

[Kd Validated] Foxk1 Rabbit Mab - P85037

Item Number: CM0011817



Introduction

Rabbit monoclonal antibody against human FOXK1, validated for WB, IP, IF/ICC, IHC-P, ChIP, ELISA. Reacts with human and mouse. Ideal for research on glucose metabolism, muscle differentiation, and autophagy.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] FOXK1 Rabbit mAb
Remarks/Alias	FOXK1L
Species	Human
GeneID (human)	221937
GeneID	221937
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 600-733 of human FOXK1 (NP_001032242.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P85037
Protein Weight	75kDa
Shipping	Ice bag