


## [Kd Validated] Gpx4 Rabbit Pab For Wb, If-P, Elisa - P36969

Item Number: CM0019726



### Introduction

Rabbit polyclonal antibody against GPX4 (glutathione peroxidase 4), validated for WB, IF-P, ELISA. Reactive with human, mouse, and rat. Suitable for ferroptosis, oxidative stress, and lipid peroxidation research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated] GPX4 Rabbit pAb                                                                                   |
| Remarks/Alias    | MCSP; SMDS; GPx-4; PHGPx; snGPx; GSHPx-4; snPHGPx; GPX4                                                          |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 2879                                                                                                             |
| GeneID           | 2879                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 30-197 of human GPX4 (NP_002076.2) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IF-P, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P36969                                                                                                           |
| Protein Weight   | 22 kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                          |