


[Kd Validated] Mrpl11 Rabbit Pab - Q9Y3B7

Item Number: CM0024605



Introduction

Rabbit polyclonal antibody targeting human MRPL11, a mitochondrial ribosomal protein. Validated for WB, IP, and ELISA. Ideal for mitochondrial protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MRPL11 Rabbit pAb
Remarks/Alias	MRPL11; CGI-113; L11MT; MRP-L11; mitochondrial ribosomal protein L11
Species	Human
GeneID (Human)	65003
GeneID	65003
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 91-192 of human MRPL11 (NP_057134.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q9Y3B7
Protein Weight	18kDa/19kDa/20kDa
Shipping	Ice bag