

[Kd Validated] Rabbit Anti-Timm13 Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Y5L4

Item Number: CM0008361

Introduction

Rabbit monoclonal antibody against human TIMM13 (Mitochondrial import inner membrane translocase subunit Tim13). Validated for WB, IF/ICC, ELISA. Suitable for mitochondrial protein import and chaperone research.

[Learn More](#)

Specification	Details
Product Name	[KD Validated] TIMM13 Rabbit mAb
Remarks/Alias	TWIST1; ACS3; BPES2; BPES3; CRS; CRS1; CSO; SCS; TWIST; bHLHa38; twist-related protein 1
Species	Human
GeneID (Human)	26517
GeneID	26517
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human TIMM13 (NP_036590.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y5L4
Protein Weight	11kDa
Shipping	Ice bag