

[Kd Validated] Stat2 Rabbit Monoclonal Antibody For Wb, Elisa - P52630

Item Number: CM0009248



Introduction

High-quality [KD Validated] STAT2 rabbit monoclonal antibody from CamelBio for WB and ELISA. Recognizes human and rat STAT2 (P52630). Ideal for interferon signaling and antiviral research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] STAT2 Rabbit mAb
Remarks/Alias	P113; IMD44; ISGF-3; PTORCH3; STAT113; [KD Validated] STAT2
Species	Human
Gene ID (Human)	6773
Gene ID	6773
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-706 of human STAT2 (NP_005410.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P52630
Protein Weight	98kDa
Shipping	Ice bag