


[Kd Validated] Trim21/Ss-A Rabbit Polyclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - P19474

Item Number: CM0019474



Introduction

Rabbit polyclonal antibody against human TRIM21/SS-A (Ro52), validated for WB, IP, IF/ICC, IHC-P, and ELISA. Recognizes human, mouse, and rat TRIM21, a 54kDa E3 ubiquitin ligase involved in innate immunity and autophagy.

[Learn More](#)

Parameter	Details
Product Name	[KD Validated] TRIM21/SS-A Rabbit pAb
Remarks/Alias	SSA; RO52; SSA1; RNF81; Ro/SSA; TRIM21/SS-A
Species	Human
GeneID (Human)	6737
GeneID	6737
Immunogen	Recombinant protein corresponding to amino acids 86-334 of human TRIM21/SS-A (NP_003132.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19474
Protein Weight	54kDa
Shipping	Ice bag