

Kdm6A Rabbit Pab - 22289

Item Number: CM0029001

Introduction

KDM6A rabbit polyclonal antibody raised against mouse KDM6A (amino acids 550-700) for Western blot, IHC-P, and ELISA; cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	KDM6A Rabbit pAb
Remarks/Alias	Utx; KDM6A
Species	Mouse
Gene ID (Human)	7403
Gene ID (Mouse)	22289
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-700 of mouse KDM6A (NP_033509.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	154kDa/157kDa
Shipping	Ice bag