

Kdm6B Rabbit Polyclonal Antibody - 216850

Item Number: CM0018690

Introduction

Kdm6b (Jmjd3) rabbit polyclonal antibody specific to mouse Kdm6b. Suitable for WB, IF/ICC, and ELISA applications. Cross-reacts with human and mouse. Protein weight 176kDa.

[Learn More](#)

Parameter	Value
Product Name	Kdm6b Rabbit pAb
Remarks/Alias	Jmjd3; 1700064E03Rik; Kdm6b
Species	Mouse
Gene ID (Human)	23135
Gene ID	216850
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-190 of mouse Kdm6b (NP_001017426.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
Protein Weight	176kDa
Shipping	Ice bag