

Kdm6B Rabbit Polyclonal Antibody - Kdm6B

Item Number: CM0015623

Introduction

Rabbit polyclonal antibody targeting mouse Kdm6b (also known as Jmjd3). Generated against a recombinant protein spanning amino acids 90-190. Suitable for Western blot and ELISA, with cross-reactivity to human and mouse Kdm6b. Predicted molecular weight: 176 kDa.

[Learn More](#)

Parameter	Value
Product Name	Kdm6b Rabbit pAb
Remarks/Alias	Jmjd3; 1700064E03Rik; Kdm6b
Species	Mouse
Gene ID (Human)	23135
Gene ID	216850
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-190 of mouse Kdm6b (NP_001017426.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
Protein Weight	176kDa
Shipping	Ice bag