

Khdrbs1/Sam68 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q07666

Item Number: CM0005880



Introduction

Rabbit monoclonal antibody targeting KHDRBS1/Sam68 (p62/p68). Validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Ideal for RNA-binding and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	KHDRBS1/Sam68 Rabbit mAb
Remarks/Alias	p62; p68; Sam68; KHDRBS1/Sam68
Species	Human
GeneID (Human)	10657
GeneID	10657
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KHDRBS1 (Q07666)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07666
Protein Weight	48kDa
Shipping	Ice bag