


Khsrp Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q92945

Item Number: CM0005696



Introduction

KHSRP Rabbit Monoclonal Antibody validated for WB, IHC-P, IP, ELISA in human, mouse, and rat. Recognizes far upstream element-binding protein 2 (FUBP2), a 73 kDa RNA-binding protein involved in mRNA splicing, transport, and decay.

[Learn More](#)

Parameter	Value
Product Name	KHSRP Rabbit mAb
Remarks/Alias	p75, FBP2, KSRP, FUBP2, KHSRP
Species	Human
GeneID (Human)	8570
GeneID	8570
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human KHSRP (Q92945)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92945
Protein Weight	73 kDa
Shipping	Ice bag