

Kiaa0513 Rabbit Polyclonal Antibody For Wb And Elisa - O60268

Item Number: CM0017876

Introduction

Rabbit polyclonal antibody targeting KIAA0513, validated for WB and ELISA in human and mouse samples. Ideal for studying cytoplasmic protein localization in brain tissues and its potential role in neuropsychiatric disorders.

[Learn More](#)

Parameter	Specification
Product Name	KIAA0513 Rabbit pAb
Remarks/Alias	KIAA0513
Species	Human
Gene ID (Human)	9764
Gene ID	9764
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-301 of human KIAA0513 (NP_001273494.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60268
Protein Weight	47kDa
Shipping	Ice bag