

Kiaa1456 Rabbit Polyclonal Antibody (Trm9L) For Wb, Elisa - Q9P272

Item Number: CM0028093

Introduction

Rabbit polyclonal antibody targeting KIAA1456 (TRM9L/tRNA methyltransferase 9B) validated for Western blot and ELISA, with cross-reactivity to human, mouse, and rat. Suitable for cancer research and tRNA modification studies.

[Learn More](#)

Parameter	Value
Product Name	KIAA1456 Rabbit pAb
Remarks/Alias	TRM9L; hTRM9L; C8orf79; KIAA1456
Species	Human
GeneID (Human)	57604
GeneID	57604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-330 of human KIAA1456 (NP_065895.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P272
Protein Weight	51kDa
Shipping	Ice bag