

Kif2B Rabbit Pab For Wb, Elisa - Q8N4N8

Item Number: CM0025761



Introduction

CamelBio's KIF2B rabbit polyclonal antibody (Q8N4N8) detects human KIF2B kinesin in WB and ELISA. Validated in human, mouse, and rat, with a predicted weight of 76 kDa. Ideal for cell division and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	KIF2B Rabbit pAb
Remarks/Alias	KIF2B
Species	Human
GeneID (human)	84643
GeneID	84643
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 544-673 of human KIF2B (NP_115948.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N4N8
Protein Weight	76kDa
Shipping	Ice bag