

Klf17 Rabbit Polyclonal Antibody For Wb And Elisa - Q5Jt82

Item Number: CM0015180

Introduction

Rabbit polyclonal antibody targeting human KLF17, validated for WB and ELISA, with cross-reactivity to mouse. Suitable for studies on transcriptional regulation, epithelial-mesenchymal transition, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	KLF17 Rabbit pAb
Remarks/Alias	ZLF393; ZNF393; Zfp393; KLF17
Species	Human
Gene ID (Human)	128209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-389 of human KLF17 (NP_775755.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5JT82
Protein Weight	43kDa
Shipping	Ice bag