

Klf6 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99612

Item Number: CM0016430

Introduction

KLF6 rabbit polyclonal antibody for detection of human KLF6 (Zf9/COPEB), a zinc-finger transcription factor. Validated in WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse. Sourced from rabbit, shipped on ice bags.

[Learn More](#)

Parameter	Value
Product Name	KLF6 Rabbit pAb
Remarks/Alias	GBF; ZF9; BCD1; CBA1; CPBP; PAC1; ST12; COPEB; KLF6
Species	Human
Gene ID (Human)	1316
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-200 of human KLF6 (NP_001291.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99612
Protein Weight	32kDa
Shipping	Ice bag