

KLHL9 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9P2J3

Item Number: CM0014287

Introduction

Rabbit polyclonal antibody raised against human KLHL9, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for mitotic progression and cytokinesis research.

[Learn More](#)

Parameter	Value
Product Name	KLHL9 Rabbit pAb
Remarks/Alias	KLHL9
Species	Human
GeneID (Human)	55958
GeneID	55958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human KLHL9 (NP_061335.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2J3
Protein Weight	69kDa
Shipping	Ice bag