

## KLK10 (Nes1) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43240

Item Number: CM0026177



### Introduction

Recombinant rabbit polyclonal antibody targeting human Kallikrein-10 (KLK10/NES1), validated for WB, IF/ICC, and ELISA. Detects human, mouse, rat KLK10. Tumor suppressor antigen for oncology research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | KLK10 Rabbit pAb                                                                                                  |
| Remarks/Alias    | NES1; PRSSL1; KLK10                                                                                               |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 5655                                                                                                              |
| GeneID           | 5655                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 43-276 of human KLK10 (NP_665895.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IF/ICC, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | O43240                                                                                                            |
| Protein Weight   | 30kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |