


Kmt2A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q03164

Item Number: CM0016140



Introduction

KMT2A rabbit polyclonal antibody validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Targets histone methyltransferase MLL1 involved in H3K4 methylation and transcriptional regulation, relevant to leukemia and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	KMT2A Rabbit pAb
Remarks/Alias	HRX; MLL; ALL1; GAS7; HTRX; MLL1; TRX1; ALL-1; CXXC7; HTRX1; MLL1A; WDSTS; KMT2A
Species	Human
GeneID (Human)	4297
GeneID	4297
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2829-2883 of human KMT2A (NP_005924.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03164
Protein Weight	432kDa
Shipping	Ice bag