

[KO Validated] Ak2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P54819

Item Number: CM0026977



Introduction

Rabbit polyclonal antibody against human AK2, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Adenylate kinase 2 is a mitochondrial enzyme critical for cellular energy homeostasis and hematopoiesis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] AK2 Rabbit pAb
Remarks/Alias	ADK2; K2
Species	Human
GeneID	204
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-232 of human AK2 (NP_001616.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54819
Protein Weight	26kDa
Shipping	Ice bag