

[KO Validated] Anti-A-Catenin Monoclonal Antibody For Wb - P35221

Item Number: CM0002574



Introduction

Rabbit monoclonal antibody against human α -Catenin (Alpha E-catenin, CTNNA1). KO validated for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Ideal for cell adhesion and junction research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] α -Catenin Rabbit mAb
Remarks/Alias	MDBS2; MDPT2; CAP102; in
Species	Human
Gene ID (Human)	1495
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human α -Catenin (P35221)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35221
Protein Weight	100kDa
Shipping	Ice bag