

Ko Validated Anti-Abcf2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ug63

Item Number: CM0017938

Introduction

KO Validated anti-ABCF2 rabbit polyclonal antibody validated in Western blot and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for ABC transporter research and diagnostic assay development, sourced from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ABCF2 Rabbit pAb
Remarks/Alias	ABC28; HUSSY18; HUSSY-18; EST133090; F2
Species	Human
Gene ID	10061 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ABCF2 (NP_005683.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UG63
Protein Weight	71kDa
Shipping	Ice bag