


[KO Validated] Anti-Aif Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - O95831

Item Number: CM0026153



Introduction

Reliable KO-validated AIF rabbit polyclonal antibody for WB, IF/ICC, IHC-P, ELISA. Detects human, mouse, rat AIF. Ideal for apoptosis and mitochondrial research. Target MW 67kDa. Swiss-Prot: O95831.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] AIF Rabbit pAb
Remarks/Alias	AIF; AUNX1; CMT2D; CMTX4; COWCK; DFNX5; NADMR; NAMS; PDCD8; COXPD6; SEMDHL; IF
Species	Human
Gene ID (Human)	9131
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 334-613 of human AIF (NP_004199.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95831
Protein Weight	67kDa
Shipping	Ice bag