

[Ko Validated] Anti-Ajuba Rabbit Polyclonal Antibody For Wb, Elisa - Q96If1

Item Number: CM0023880

Introduction

Rabbit polyclonal antibody against human Ajuba, validated by knockout. Suitable for WB and ELISA; cross-reacts with human and rat. For research on cell adhesion, signal transduction, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Ajuba Rabbit pAb
Remarks/Alias	JUB; ba
Species	Human
Gene ID	84962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 201-290 of human Ajuba (NP_116265.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96IF1
Protein Weight	57kDa
Shipping	Ice bag