

[KO Validated] Anti-Ak2 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P54819

Item Number: CM0002691



Introduction

[KO Validated] Recombinant rabbit monoclonal antibody targeting human adenylate kinase 2 (AK2). Suitable for Western blot, immunofluorescence/ICC, and ELISA. Recognizes human, mouse, and rat AK2. AK2 is essential for cellular energy homeostasis and hematopoiesis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] AK2 Rabbit mAb
Remarks/Alias	ADK2; K2
Species	Human
GeneID (Human)	204
GeneID	204
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human AK2 (P54819).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54819
Protein Weight	26 kDa
Shipping	Ice bag