

[KO Validated] Anti-Alpha Catenin Rabbit Polyclonal Antibody - P35221

Item Number: CM0029706



Introduction

KO validated rabbit polyclonal antibody against human alpha-catenin (CTNNA1) for WB, IHC-P, IP, ELISA; cross-reacts with mouse and rat; ideal for cell adhesion and junction research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] α -Catenin Rabbit pAb
Remarks/Alias	MDBS2; MDPT2; CAP102; in
Species	Human
GeneID (Human)	1495
GeneID	1495
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 637-906 of human α -Catenin (NP_001894.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35221
Protein Weight	100kDa
Shipping	Ice bag