


Ko Validated Anti-Alpha-Synuclein Monoclonal Antibody For Wb, Ihc, If, Elisa - P37840

Item Number: CM0011842



Introduction

Mouse monoclonal antibody against human Alpha-synuclein, validated for Western blot, IHC, IF, and ELISA. Cross-reacts with mouse and rat. Ideal for Parkinson's disease and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Alpha-synuclein Mouse mAb
Remarks/Alias	PD1; NACP; PARK1; PARK4
Species	Human
Gene ID	6622
Immunogen	Recombinant fusion protein corresponding to amino acids 1-140 of human Alpha-synuclein (NP_000336.1)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37840
Protein Weight	14kDa
Shipping	Ice bag