

[KO Validated] Anti-Ampk α 1 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - Q13131

Item Number: CM0008779



Introduction

KO-validated rabbit monoclonal antibody targeting human AMPK α 1 (Gene ID: 5562). Validated for WB, IP, IF/ICC, IHC-P, ELISA applications and cross-reacts with human, mouse, rat. Essential for studying cellular energy metabolism.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] AMPK α 1 Rabbit mAb
Remarks/Alias	AMPK; AMPK α 1; AMPK alpha 1
Species	Human
GeneID (Human)	5562
GeneID	5562
Immunogen	Synthetic peptide corresponding to amino acids 381-480 of human AMPK α 1 (NP_006242.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13131
Protein Weight	64kDa
Shipping	Ice bag