


[KO Validated] Anti-AMPK α 2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P54646

Item Number: CM0028623



Introduction

KO Validated Anti-AMPK α 2 rabbit polyclonal antibody (pAb) for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Detects 62 kDa AMPK alpha-2 catalytic subunit. Ideal for energy metabolism and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] AMPK α 2 Rabbit pAb
Remarks/Alias	AMPK; AMPK2; PRKAA; AMPK α 2; α 2
Species	Human
GeneID (Human)	5563
GeneID	5563
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54646
Protein Weight	62kDa
Shipping	Ice bag