

Ko Validated Anti-Ankrd52 Polyclonal Antibody For Wb, Elisa - Q8Nb46

Item Number: CM0021164

Introduction

KO validated rabbit polyclonal antibody targeting ANKRD52, a regulatory subunit of protein phosphatase 6. Suitable for WB and ELISA, and cross-reacts with human, mouse, and rat ANKRD52.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ANKRD52 Rabbit pAb
Remarks/Alias	ANKRD33; 52
Species	Human
GeneID	283373
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-660 of human ANKRD52 (NP_775866.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NB46
Protein Weight	115kDa
Shipping	Ice bag