


[KO Validated] Anti-Arhgef1 Rabbit Polyclonal Antibody - Q92888

Item Number: CM0026215



Introduction

KO validated rabbit polyclonal antibody targeting human ARHGEF1, a Rho guanine nucleotide exchange factor. Suitable for WB, IHC-P, ELISA. Detects endogenous protein. Cross-reacts with human and rat.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ARHGEF1 Rabbit pAb
Remarks/Alias	LSC; GEF1; IMD62; LBCL2; SUB1.5; P115-RHOGEF; F1
Species	Human
GeneID (Human)	9138
GeneID	9138
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 783-912 of human ARHGEF1 (NP_004697.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q92888
Protein Weight	102kDa
Shipping	Ice bag