


[KO Validated] Anti-Atg4A Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wyn0

Item Number: CM0024292



Introduction

Rabbit polyclonal antibody targeting ATG4A, a cysteine protease essential for autophagy. Validated by knockout, suitable for Western blot and ELISA applications. Reacts with human and mouse. Targets the C-terminal region (aa 289-398).

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ATG4A Rabbit pAb
Remarks/Alias	APG4A; AUTL2; HsAPG4A; 4A
Species	Human
GeneID (Human)	115201
GeneID	115201
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 289-398 of human ATG4A (NP_443168.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WYN0
Protein Weight	45kDa
Shipping	Ice bag