

Ko Validated Anti-Atg4D Rabbit Polyclonal Antibody For Wb And Elisa - Q86TL0

Item Number: CM0020228

Introduction

KO-validated anti-ATG4D rabbit polyclonal antibody. Detects a 53 kDa band in WB; also suitable for ELISA. Recognizes human, mouse, and rat ATG4D. Research use in autophagy and mitochondrial studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ATG4D Rabbit pAb
Remarks/Alias	APG4D; AUTL4; APG4-D; HsAPG4D; 4D
Species	Human
Gene ID	84971
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human ATG4D (NP_116274.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86TL0
Protein Weight	53kDa
Shipping	Ice bag