


Ko Validated Anti-Bad Rabbit Polyclonal Antibody For Wb And Elisa - Q92934

Item Number: CM0020486



Introduction

Rabbit polyclonal antibody against human Bad (Bcl2-associated agonist of cell death), validated in WB and ELISA. Cross-reacts with mouse. Targets the pro-apoptotic Bcl-2 family member involved in mitochondrial-mediated cell death.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Bad Rabbit pAb
Remarks/Alias	BBC2; BCL2L8; ad
Species	Human
GeneID (human)	572
GeneID	572
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 68-168 of human Bad (NP_004313.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92934
Protein Weight	18kDa
Shipping	Ice bag