


[KO Validated] Anti-Bax Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q07812

Item Number: CM0002404



Introduction

Rabbit monoclonal antibody against human Bax, KO validated, for WB, IHC-P, ELISA. Recognizes human, mouse, and rat Bax. Useful for apoptosis research. Immunogen: synthetic peptide from human Bax (Q07812).

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Bax Rabbit mAb
Remarks/Alias	BCL2L4; Bax
Species	Human
GenelD (Human)	581
GenelD	581
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07812
Protein Weight	21kDa
Shipping	Ice bag