


KO Validated Anti-Bax Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q07812

Item Number: CM0002312



Introduction

KO validated rabbit monoclonal antibody against human Bax, suitable for WB, IHC-P, IF/ICC, IP, and ELISA. Key reagent for apoptosis research, targeting the mitochondrial outer membrane protein Bax.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Bax Rabbit mAb
Remarks/Alias	BCL2L4; ax
Species	Human
Gene ID	581
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 10-70 of human Bax (NP_620116.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	Q07812
Protein Weight	21kDa
Shipping	Ice bag