


[KO Validated] Anti-Brd7 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Npi1

Item Number: CM0020805



Introduction

KO validated rabbit polyclonal antibody targeting BRD7 for WB and ELISA. Validated in human, BRD7 is a bromodomain-containing protein involved in chromatin remodeling and Wnt signaling.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] BRD7 Rabbit pAb
Remarks/Alias	BP75; NAG4; CELTIX1; SMARCI1; D7
Species	Human
GeneID (Human)	29117
GeneID	29117
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human BRD7 (NP_037395.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NPI1
Protein Weight	74kDa
Shipping	Ice bag