

KO Validated Anti-Brg1/Smarca4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P51532

Item Number: CM0020613



Introduction

KO validated rabbit polyclonal antibody against BRG1/SMARCA4, suitable for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with human, mouse, rat. For chromatin remodeling and transcription research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] BRG1/SMARCA4 Rabbit pAb
Remarks/Alias	BRG1; CSS4; SNF2; SWI2; MRD16; RTPS2; BAF190; SNF2L4; SNF2LB; hSNF2b; BAF190A; SNF2-beta; A4
Species	Human
GeneID (Human)	6597
Gene ID	6597
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 680-770 of human BRG1/SMARCA4 (NP_003063.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51532
Protein Weight	185kDa
Shipping	Ice bag