

[Ko Validated] Anti-Calnexin Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P27824

Item Number: CM0012533

Introduction

KO-validated rabbit monoclonal antibody targeting human Calnexin (P27824), an ER chaperone involved in glycoprotein folding and quality control. Validated for WB, IF/ICC, IHC-P, and ELISA.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Calnexin Rabbit mAb
Remarks/Alias	CNX; P90; IP90
Species	Human
GeneID (Human)	821
GeneID	821
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Calnexin (NP_001737.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P27824
Protein Weight	68kDa
Shipping	Ice bag