


[KO Validated] Anti-Caspase-8 Rabbit Monoclonal Antibody For Wb, Elisa - Q14790

Item Number: CM0003845



Introduction

KO validated rabbit monoclonal antibody targeting human Caspase-8 (Q14790). Suitable for Western blot and ELISA applications, it specifically detects the ~55 kDa apoptosis, necroptosis, and pyroptosis initiator protease, supporting IVD research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Caspase-8 Rabbit mAb
Remarks/Alias	CAP4; MACH; MCH5; FLICE; ALPS2B; Casp-8; Caspase-8
Species	Human
Gene ID (Human)	841
Gene ID	841
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Caspase-8 (Q14790).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14790
Protein Weight	55kDa
Shipping	Ice bag