


Ko Validated Anti-Caspase-9 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P55211

Item Number: CM0003394



Introduction

KO Validated anti-Caspase-9 rabbit monoclonal antibody for WB, IF/ICC, ELISA in human, mouse, rat. Detects 46kDa protein; targets apoptotic protease Mch-6.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Caspase-9 Rabbit mAb
Remarks/Alias	MCH6; APAF3; APAF-3; PPP1R56; ICE-LAP6; Caspase-9
Species	Human
GeneID (Human)	842
GeneID	842
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 139-416 of Caspase-9 (NP_001220.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55211
Protein Weight	46kDa
Shipping	Ice bag