

[Ko Validated] Anti-Cd133 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - O43490

Item Number: CM0015973

Introduction

Validated rabbit polyclonal antibody against human CD133 (Prominin-1). Suitable for WB, IF/ICC, IP, and ELISA applications; cross-reacts with mouse and rat. Ideal for stem cell and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CD133 Rabbit pAb
Remarks/Alias	RP41; AC133; CD133; MCDR2; STGD4; CORD12; PROML1; MSTP061; 33
Species	Human
GeneID (Human)	8842
GeneID	8842
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43490
Protein Weight	97kDa
Shipping	Ice bag