

Ko Validated Anti-Cd147/Bsg Polyclonal Antibody For Wb, Elisa - P35613

Item Number: CM0019524



Introduction

KO validated rabbit polyclonal antibody against human CD147/BSG (Basigin, EMMPRIN). Suitable for WB and ELISA, targeting a transmembrane glycoprotein involved in tumor invasion, angiogenesis, and microbial infection.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CD147/BSG Rabbit pAb
Remarks/Alias	OK; 5F7; TCSF; CD147; EMPRIN; HAB18G; EMMPRIN; SG
Species	Human
Gene ID (Human)	682
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-200 of human CD147/BSG (NP_940991.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P35613
Protein Weight	42kDa
Shipping	Ice bag