

[Ko Validated] Anti-Cd47 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q61735

Item Number: CM0005906

Introduction

KO validated rabbit monoclonal antibody targeting mouse CD47 (IAP), ideal for WB, IHC-P, ELISA. Recognizes leukocyte surface antigen involved in immune regulation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CD47 Rabbit mAb
Remarks/Alias	IAP; Itgp; 9130415E20Rik; B430305P08Rik
Species	Mouse
Gene ID (Human)	961
Gene ID	16423
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-140 of mouse CD47 (NP_034711.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	Q61735
Protein Weight	33 kDa/35 kDa
Shipping	Ice bag