

[KO Validated] Anti-Cdk4 Rabbit Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - P11802

Item Number: CM0002256



Introduction

KO validated rabbit monoclonal antibody targeting human CDK4 (P11802), suitable for WB, IP, IHC-P, ELISA. Cross-reacts with human, mouse, rat, monkey. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDK4 Rabbit mAb
Remarks/Alias	CMM3; PSK-J3; K4
Species	Human
GeneID (Human)	1019
GeneID	1019
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P11802
Protein Weight	34kDa
Shipping	Ice bag