


[KO Validated] Anti-Cdk8 Polyclonal Antibody For Wb, Ip, Elisa - P49336

Item Number: CM0024776



Introduction

CamelBio supplies [KO Validated] CDK8 Rabbit Polyclonal Antibody (P49336) validated for WB, IP, ELISA. Reliable for transcription regulation research in human, mouse, rat samples.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDK8 Rabbit pAb
Remarks/Alias	K35, IDDHBA, K8
Species	Human
Gene ID	1024
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 365-464 of human CDK8 (NP_001251.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49336
Protein Weight	53kDa
Shipping	Ice bag