

[KO Validated] Anti-Cdkn1A/P21 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P38936

Item Number: CM0003039



Introduction

KO validated rabbit monoclonal antibody targeting human CDKN1A/p21 (P38936). Suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for cell cycle, DNA damage, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDKN1A/p21 Rabbit mAb
Remarks/Alias	P21; CIP1; SDI1; WAF1; CAP20; CDKN1; MDA-6; p21CIP1; 21
Species	Human
GeneID (Human)	1026
GeneID	1026
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-164 of human CDKN1A/p21CIP1 (NP_000380.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38936
Protein Weight	18kDa
Shipping	Ice bag