

[KO Validated] Anti-Cdkn1A/P21 Rabbit Polymab® Monoclonal Antibody For Wb, If/Icc, Elisa - P38936 / P39689

Item Number: CM0006049



Introduction

Rabbit recombinant monoclonal antibody targeting CDKN1A/p21, validated for WB, IF/ICC, ELISA. Reacts with human and mouse, 18 kDa. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDKN1A/p21 Rabbit PolymAb®
Remarks/Alias	P21; CIP1; SDI1; WAF1; CAP20; CDKN1; MDA-6; p21CIP1
Species	Human/Mouse
Gene ID (Human)	1026
Gene ID	1026/12575
Immunogen	Recombinant protein of human/mouse CDKN1A/p21
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P38936/P39689
Protein Weight	18 kDa
Shipping	Ice bag